KIT LOT NUMBER: B-1662-2
EXPIRY DATE: 2014/02

REFERENCE PLASMA VALUE: **97.5%**

ABNORMAL CONTROL RANGE:
(21.6%- 32.4%)

abp

Ristocetin Co-Factor Assay Kit

Ristocetin-Kofaktor-Testkit

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
Ristocetin Co-Factor Assay Kit

INTENDED PURPOSE

The abp Ristocetin co-factor assay kit is intended for the quantitation of von Willebrand factor activity through the reflective assay of ristocetin cofactor activity.

TEST PRINCIPLE

The abp Ristocetin co-factor assay kit measures the ability of a patients plasma to agglutinate formalin-fixed platelets in the presence of ristocetin. Serial dilutions of a reference plasma with an assigned ristocetin co-factor value is performed. A standard curve is devised. Patient plasmas are assayed in the presence of formalin-fixed platelets and ristocetin to determine %/min agglutination. Results can be interpolated against standard curve to determine ristocetin co-factor activity.

MATERIALS REQUIRED BUT NOT SUPPLIED

Pipettes for 500ul, 200ul and 25ul.
Aggregometer with chart recorder or plotting software
Aggregometer cuvettes
Plastic tubes for dilution of standard and patients

PRESENTATION

abp Ristocetin co-factor assay kit is supplied as:
4 x 5ml lyophilised formalin-fixed human platelets
2 x 1ml 10mg/ml Ristocetin lyophilised
2 x 1.0ml Reference standard derived from normal human plasma - lyophilised
2 x 0.5ml Abnormal control derived from human plasma
2 x 25ml Tris-buffered saline with preservatives and stabilisers

REAGENT PREPARATION

Lyophilised Platelets

Allow vials to equilibrate to room temperature before reconstitution. Reconstitute with 5ml of Tris-buffered saline. Leave to stand at least 15 minutes. Mix gently before use by inverting vial. Store at 2–8°C. Stable for 56 days upon reconstitution at 2 - 8°C. Stable for 7 hours upon reconstitution at room temperature.

Ristocetin

Allow vial to equilibrate to room temperature before use. Reconstitute with 1ml of reagent grade water to achieve a final concentration of 10mg/ml. Allow to stand for 15 minutes. Ensure all particulate matter has dissolved before use. Store at 2-8°C. Any unused portion can be frozen and stored at -20°C. It is also stable –20°C for 14 days.

Reference Standard

Allow vial to equilibrate to room temperature before reconstitution. Reconstitute with 1.0ml of reagent grade water. Allow to stand for 10–15 minutes with occasional swirling. Ensure all particulate matter is dissolved before use. Store at 2-8°C. Stable until expiry on vial label. Reconstituted plasma is stable for 7 hours at room temperature

Abnormal control patient plasma

Allow vial to equilibrate to room temperature before reconstitution. Reconstitute with 0.5ml of reagent grade water. Allow to stand for 10–15 minutes with occasional swirling. Ensure all particulate matter is dissolved before use. Store at 2-8°C. Stable until expiry on vial label. Reconstituted plasma is stable for 7 hours at room temperature.

Tris-Buffered Saline

Buffer is ready prepared. Contains sodium azide 0.1% and Thiomersal 0.02%. Store at 2-8°C. Stable until expiry on vial label.

WARNINGS AND PRECAUTIONS

For in-vitro diagnostic use only.

Each unit of plasma used in this product has been screened for the presence of antibody to Human Immunodeficiency Virus (HIV) Type I and II, Hepatitis B surface antigen and Hepatitis C and found negative.

However no test can offer complete assurance that human plasma will not transmit infectious disease. As with all materials of human origin, this product should be treated as infectious and necessary precautions should be implemented. Dispose of waste material according to individual countries regulations.

SPECIMEN PREPARATION

Specimen: Plasma obtained from whole blood venous collection and 3.2% sodium citrate as an anticoagulant is the specimen of choice.

Specimen collection and preparation: Obtain venous blood sample by venipuncture. Immediately transfer the blood into a tube containing suitable anti-coagulant (Sodium citrate 3.2%) and mix well by inversion. Centrifuge specimen at 1500g for 20 minutes. Remove plasma from tubes within 60 minutes. Test plasma sample within 4 hours for optimum results or store frozen at -20°C for up to 2 weeks. Avoid freeze thawing technique.(NCCLS Standard H21-A2 ')

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TEST PROCEDURE

Assay calibration

Prepare the standard curve by making serial dilutions of reference plasma in buffer as follows:

Tube	Vol of Buffer	Vol of Reference	%Assigned Reference Value
1	0.3 ml	0.3ml	100
2	0.1 ml	0.3 ml Tube 1	75
3	0.1 ml	0.2 ml Tube 2	50
4	0.2 ml	0.2 ml Tube 3	25

Prepare test plasma dilution by adding 0.1ml Test plasma to 0.1ml TBS in a plastic tube. Mix well.

To prepare a blank needed to set the 100% baseline for each sample cuvette, add 0.10ml of TBS to 0.05ml of reconstituted platelets.

Assay procedure

1. Add 0.2ml reconstituted platelets to an aggregometer cuvette, ensuring stir bar is present.
2. Add 0.025ml Ristocetin to the sample cuvette. Mix well. Incubate for 120 seconds.
3. During incubation, set the 100% baseline by using the blank made in section previously. If a 0% baseline is required then TBS is normally used.
4. Begin aggregation by adding 0.025ml 100% reference standard dilution to Ristocetin containing platelets. Observe and record aggregation pattern, ending agglutination at 180 seconds.
5. Repeat steps 1 – 4 for all reference standard dilutions 75%, 50% and 25%.
6. Repeat steps 1 - 4 for patient test sample (1+1 dilution).
7. Abnormal samples (below 40%) should be repeated using undiluted plasma and divided by 2 to obtain higher degree of result precision.

QUALITY CONTROL

Absorbed von Willebrand control – abnormal is included in the test kit. The Ristocetin Cofactor value of this control has a set range. When assayed as a test plasma and interpolated against assay standard curve, the given result should fall within this range. If the abnormal control does not fall into range, repeat standard curve and control. If control does not fall into range again, this is indicative of reagent deterioration. Contact abp for further assistance.

RESULTS

A. Determination of slope

1. See instrument manufacturer's instructions.

B. Preparation of standard curve from slope values

1. Using log-log graph paper supplied, plot the % reference plasma 100%, 75%, 50% and 25% adjusting actual % to reference plasma ristocetin co-factor activity % level. Plot ristocetin % on x axis and slope value on the y axis.
2. Draw a line of best fit through the four points.

C. Determine Ristocetin Co-factor activity in patient samples

1. From the slope value of the patient plasma, interpolate against line of best fit, extending line to x axis and establishing the level of ristocetin Co-factor activity of the patient. It is also possible to use automated plotting software to determine Ristocetin Co-factor activity levels of patient samples.

EXPECTED VALUES

In general three types of vWd are recognised. In general type I vWd is characterised by reduced levels of all the normal multimers in plasma and platelets. Type II shows an abnormal pattern of plasma multimers with normal or abnormal platelet patterns, and type III have undetectable vWF. Normal results vary between 63% and 99%. In von Willebrands disease results are generally below 40% and usually below 30%².

LIMITATIONS

Situations such as Pregnancy, infusion of commercial Factor VIII concentrates³ may fail to reflect accurately von Willebrand's disease. However Ristocetin Co-factor activity has become perhaps the most reliable indicator of von Willebrands disease. It appears to be the assay most closely related to the biological activity of von Willebrand factor.

RISK AND SAFETY

Harmful contains sodium azide.

Normal laboratory health and safety procedures should be adhered to.

REFERENCES

1. NCCLS, Collection, transport and processing of blood specimens for coagulation testing and performance of coagulation assays, 2nd edition H21-A3 1991
2. Thompson, J.M., Blood coagulation and haemostasis, a practical guide, third edition, Longman group (FE) pg. 142,145,192, 1985.
3. Blatt, P.M., Brinkhous, K.M.M., Culp, H.R., *et al.*, Antithrombotic factor concentrate therapy in von Willebrands disease, J Am Med Assn, 236:2770-2772, 1976.

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RÉFÉRENCES

1. NCCLS, Collection, transport and processing of blood specimens for coagulation testing and performance of coagulation assays, 2nd edition H21-A3 1991
2. Thompson, J.M., Blood coagulation and haemostasis, a practical guide, third edition, Longman group (FE) pg. 142,145,192, 1985.
3. Blatt, P.M., Brinkhous, K.M.M., Culp, H.R., *et al.*, Antithrombotic factor concentrate therapy in von Willebrand's disease, J Am Med Assn, 236:2770-2772, 1976.